Desktop LIBS analyzer

This device can be placed on the table and next to the production area. Then we easily open the door of the device and place the sample on the sample holder and move it to the desired location using the laser level. Finally, we close the door and the measurement begins.

This device is used to identify elements in materials such as coal, steel slag, gems, food and medicine, etc.

Comparing LIBS to other common methods

Method	LIBS	OES	XRF
Destructive effect	Yes (small)	Yes	No
Elements	Lithium, Beryllium, Magnesium, Aluminum, Silicon	Carbon, Boron, Phosphorus, Nitrogen, Sulfur	Tungsten (heavy metal)
Test time	1-3 s	Laser: 6 s Spark: The slowest technique	5-20 s
Harmful radiation	No	No	Yes
Maintenance and repair conditions	Easy maintenance, low repair cost, durable	Laser: Low maintenance Spark: Regular cleaning and expensive repairs	Expensive repairs

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Technical specifications of the device

- Easy to use
- Quick detection
- On the spot detection
- High accuracy of detection
- Very little destructive effect
- Need less amount of sample
- No need for sample preparation
- Detection of all available light and heavy elements
- Possibility of customization for one or more specific elements

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
Technical specifications of the device

Device dimensions	20x23x19 inches
Device weight	45 kg
Limit of detection (LOD)	10 ppm
Spectrometer channels	4 to 7
Detection time	20 s
Maximum sample size	5x9x5 inches
Spot size	20-200 μm
Laser penetration depth	1-100 μm
Laser type	1064 nm
Laser energy	Up to 200 mJ/pulse
Software	PartoSoft

Physical specifications of the device

- Detection of elements' spectral lines
- Calibrated with mercury lamp
- 120 volts AC and 10 amp
- The radius of the destruction spot is 50 μm to 2

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